

Product Datasheet

TIMD4 antibody (orb166253)

Description	Rabbit polyclonal antibody to TIMD4 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TIMD4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	41578
Uniprot ID	Q96H15
NCBI	NP_001140198.1 , NP_612388.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

