

Product Datasheet

Mouse Lysyl oxidase homolog 3 ELISA Kit (orb1662624)

Description	Mouse Lysyl oxidase homolog 3 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78-5000 pg/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 T. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.