

Product Datasheet

Mouse Kinesin Family, Member 5B (KIF5B) ELISA Kit (orb1662625)

Description	Mouse Kinesin Family, Member 5B (KIF5B) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78.1-5000 pg/mL. Sensitivity: 29.4 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates or other biological fluids.
Expiration Date	Please enquire.