

Product Datasheet

Mouse V-Jun Sarcoma Virus 17 Oncogene Homolog (JUN) ELISA Kit (orb1662629)

Description	Mouse V-Jun Sarcoma Virus 17 Oncogene Homolog (JUN) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.391-25 ng/mL. Sensitivity: 0.172 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.