

Product Datasheet

Human V-Jun Sarcoma Virus 17 Oncogene Homolog (JUN) ELISA Kit (orb1662630)

Description Human V-Jun Sarcoma Virus 17 Oncogene Homolog (JUN) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 1.56-100 ng/mL. Sensitivity: 0.61 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.

Expiration Date Please enquire.