

## Product Datasheet

### Mouse JNK1/MAPK8 Associated Membrane Protein (JKAMP) ELISA Kit (orb1662631)

|                            |                                                                                                                                                                                              |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse JNK1/MAPK8 Associated Membrane Protein (JKAMP) ELISA Kit                                                                                                                               |
| <b>Reactivity</b>          | Mouse                                                                                                                                                                                        |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                        |
| <b>Note</b>                | For research use only                                                                                                                                                                        |
| <b>Application notes</b>   | Detect Range: 0.781-50 ng/mL. Sensitivity: 0.26 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates or other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                              |