

Product Datasheet

ST6GALNAC5 antibody (orb166266)

Description	Rabbit polyclonal antibody to ST6GALNAC5 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ST6GALNAC5
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	38443
Uniprot ID	Q9BVH7
NCBI	NP_112227.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

