

Product Datasheet

Mouse Cationic Amino Acid Transporter 2 (CAT2) ELISA Kit (orb1662685)

Description	Mouse Cationic Amino Acid Transporter 2 (CAT2) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.124 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates or other biological fluids.
Expiration Date	Please enquire.