

Product Datasheet

Mouse ATPase, Na⁺/K⁺ Transporting Beta 4 Polypeptide (ATP1b4) ELISA Kit (orb1662696)

Description	Mouse ATPase, Na ⁺ /K ⁺ Transporting Beta 4 Polypeptide (ATP1b4) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78.1-5000 pg/mL. Sensitivity: 27.8 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates or other biological fluids.
Expiration Date	Please enquire.