

Product Datasheet

Mouse Aquaporin 2, Collecting Duct (AQP2) ELISA Kit (orb1662714)

Description	Mouse Aquaporin 2, Collecting Duct (AQP2) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.059 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.