

Product Datasheet

Mouse Amyloid Beta Precursor Protein Binding Protein B3 (APBB3) ELISA Kit (orb1662725)

Description	Mouse Amyloid Beta Precursor Protein Binding Protein B3 (APBB3) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.625-40 ng/mL. Sensitivity: 0.255 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates or other biological fluids.
Expiration Date	Please enquire.