

Product Datasheet

Human Aldehyde Dehydrogenase 7 Family, Member A1 (ALDH7A1) ELISA Kit (orb1662741)

Description Human Aldehyde
Dehydrogenase 7 Family,
Member A1 (ALDH7A1) ELISA Kit

Reactivity Human

**Tested
Applications** ELISA

Note For research use only

Application notes Detect Range: 0.781-50 ng/mL.
Sensitivity: 0.29 ng/mL. Assay
Time: 1-4.5 h. Detection
wavelength: 450 nm. Sample
Volume: 100 µL. Sample Type:
serum, plasma, tissue
homogenates, cell lysates, cell
culture supernates or other
biological fluids.

Expiration Date Please enquire.