

Product Datasheet

Rat Alcohol Dehydrogenase 7 (ADH7) ELISA Kit (orb1662747)

Description Rat Alcohol Dehydrogenase 7 (ADH7)
ELISA Kit

Reactivity Rat

Tested ELISA

Applications

Note For research use only

Application notes Detect Range: 0.156-10 ng/mL.
Sensitivity: 0.058 ng/mL. Assay
Time: 1-4.5 h. Detection wavelength:
450 nm. Sample Volume: 100 µL.
Sample Type: serum, plasma, tissue
homogenates or other biological
fluids.

Expiration Date Please enquire.