

Product Datasheet

ZNF205 antibody (orb166275)

Description	Rabbit polyclonal antibody to ZNF205 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ZNF205
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	60630
Uniprot ID	O95201
NCBI	NP_003447.2 , NP_001265087.1 , NP_001035893.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

