

Product Datasheet

Rat Zona Pellucida Glycoprotein 2, Sperm Receptor (ZP2) ELISA Kit (orb1662754)

Description Rat Zona Pellucida Glycoprotein 2, Sperm Receptor (ZP2) ELISA Kit

Reactivity Rat

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.312-20 ng/mL. Sensitivity: 0.125 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.

Expiration Date Please enquire.