

Product Datasheet

Human Chemokine C-Motif Ligand 2 (XCL2) ELISA Kit (orb1662784)

Description	Human Chemokine C-Motif Ligand 2 (XCL2) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.058 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.