

Product Datasheet

Human WAP Four Disulfide Core Domain Protein 5 (WFDC5) ELISA Kit (orb1662818)

Description Human WAP Four Disulfide Core Domain Protein 5 (WFDC5) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 78.1-5000 pg/mL. Sensitivity: 35 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.

Expiration Date Please enquire.