

Product Datasheet

Mouse Wiskott Aldrich Syndrome Protein Family, Member 2 (WASF2) ELISA Kit (orb1662820)

Description Mouse Wiskott Aldrich Syndrome Protein Family, Member 2 (WASF2) ELISA Kit

Reactivity Mouse

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.156-10 ng/mL. Sensitivity: 0.053 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.

Expiration Date Please enquire.