

Product Datasheet

Human Von Willebrand Factor A Domain Containing Protein 3A (vWA3A) ELISA Kit (orb1662826)

Description Human Von Willebrand Factor A
Domain Containing Protein 3A
(vWA3A) ELISA Kit

Reactivity Human

**Tested
Applications** ELISA

Note For research use only

Application notes Detect Range: 0.468-30 ng/mL.
Sensitivity: 0.15 ng/mL. Assay
Time: 1-4.5 h. Detection
wavelength: 450 nm. Sample
Volume: 100 µL. Sample Type:
plasma.

Expiration Date Please enquire.