

Product Datasheet

OVCH1 antibody (orb166292)

Description	Rabbit polyclonal antibody to OVCH1 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	OVCH1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	125066
Uniprot ID	Q7RTY7
NCBI	NP_899234.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

