

Product Datasheet

LRRC8A antibody (orb166293)

Description	Rabbit polyclonal antibody to LRRC8A (C-term).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	LRRC8A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	94199
Uniprot ID	Q8IWT6
NCBI	NP_062540.2 , NP_001120717.1 , NP_001120716.1
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.

