

Product Datasheet

Human UDP-Glucose Glycoprotein Glucosyltransferase 1 (UGGT1) ELISA Kit (orb1662937)

Description Human UDP-Glucose Glycoprotein
Glucosyltransferase 1 (UGGT1)
ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.312-20 ng/mL.
Sensitivity: 0.113 ng/mL. Assay
Time: 1-4.5 h. Detection
wavelength: 450 nm. Sample
Volume: 100 µL. Sample Type:
tissue homogenates, cell lysates
or other biological fluids.

Expiration Date Please enquire.