

Product Datasheet

Guinea pig Ubiquitin (Ub) ELISA Kit (orb1662967)

Description	Guinea pig Ubiquitin (Ub) ELISA Kit
Reactivity	Guinea pig
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 1.56-100 ng/mL. Sensitivity: 0.36 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.