

Product Datasheet

Human Ubiquitin Conjugating Enzyme E2L3 (UBE2L3) ELISA Kit (orb1662971)

Description	Human Ubiquitin Conjugating Enzyme E2L3 (UBE2L3) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.072 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates or other biological fluids.
Expiration Date	Please enquire.