

Product Datasheet

RBM7 antibody (orb166298)

Description	Rabbit polyclonal antibody to RBM7 (C-term).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	RBM7
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	30504
Uniprot ID	Q9Y580
NCBI	NP_057174.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

