

Product Datasheet

Mouse High mobility group protein B1 ELISA Kit (orb1662986)

Description	Mouse High mobility group protein B1 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 100 T. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.