

Product Datasheet

Human Transient Receptor Potential Cation Channel Subfamily M, Member 1 (TRPM1) ELISA Kit (orb1663012)

Description Human Transient Receptor Potential Cation Channel Subfamily M, Member 1 (TRPM1) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.156-10 ng/mL. Sensitivity: 0.051 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates or other biological fluids.

Expiration Date Please enquire.