

Product Datasheet

MS4A13 antibody (orb166321)

Description	Rabbit polyclonal antibody to MS4A13 (N-term).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	MS4A13
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	17307
Uniprot ID	Q5J8X5
NCBI	NP_001265249.1 , NP_001012417.2 , NP_001094379.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

