

Product Datasheet

Hras1 antibody (orb166323)

Description	Rabbit polyclonal antibody to Hras1
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	Hras
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	21298
Uniprot ID	Q61411
NCBI	NP_001123915.1 , NP_001123916.1 , NP_032310.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

