

Product Datasheet

Bovine TEK Tyrosine Kinase, Endothelial (Tie2) ELISA Kit (orb1663232)

Description	Bovine TEK Tyrosine Kinase, Endothelial (Tie2) ELISA Kit
Reactivity	Bovine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 6.25-400 ng/mL. Sensitivity: 1.56 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.