

Product Datasheet

DPCR1 antibody (orb166324)

Description	Rabbit polyclonal antibody to DPCR1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	MUCL3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	56347
Uniprot ID	Q3MIW9
NCBI	NP_543146.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

