

Product Datasheet

Rabbit Transforming Growth Factor Beta Induced Protein (TGF β I) ELISA Kit (orb1663269)

Description	Rabbit Transforming Growth Factor Beta Induced Protein (TGF β I) ELISA Kit
Reactivity	Rabbit
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.036 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.