

Product Datasheet

Monkey Transforming Growth Factor Beta 1 (TGFB1) ELISA Kit (orb1663287)

Description Monkey Transforming Growth Factor Beta 1 (TGFB1) ELISA Kit

Reactivity Monkey

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 31.2-2000 pg/mL. Sensitivity: 13.4 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell culture supernates or other biological fluids.

Expiration Date Please enquire.