

Product Datasheet

Rabbit Transforming Growth Factor Beta 1 (TGFb1) ELISA Kit (orb1663288)

Description	Rabbit Transforming Growth Factor Beta 1 (TGFb1) ELISA Kit
Reactivity	Rabbit
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6-1000 pg/mL. Sensitivity: 6.1 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, platelet-poor plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.