

Product Datasheet

Rat Taurine Transporter (TAUT) ELISA Kit (orb1663352)

Description Rat Taurine Transporter (TAUT) ELISA Kit

Reactivity Rat

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.156-10 ng/mL.
Sensitivity: 0.052 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates or other biological fluids.

Expiration Date Please enquire.