

Product Datasheet

RNF182 antibody (orb166342)

Description	Rabbit polyclonal antibody to RNF182 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	RNF182
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	27402
Uniprot ID	Q8N6D2
NCBI	NP_001158504.1 , NP_001158506.1 , NP_001158505.1 , NP_689950.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

