

Product Datasheet

ZSCAN21 antibody (orb166355)

Description	Rabbit polyclonal antibody to ZSCAN21
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ZSCAN21
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	53658
Uniprot ID	Q9Y5A6
NCBI	NP_666019.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

