

Product Datasheet

Rabbit Selectin, Endothelium (SELE) ELISA Kit (orb1663634)

Description	Rabbit Selectin, Endothelium (SELE) ELISA Kit
Reactivity	Rabbit
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 6.25-400 ng/mL. Sensitivity: 2.27 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.