

Product Datasheet

ARMCX6 antibody (orb166368)

Description	Rabbit polyclonal antibody to ARMCX6 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ARMCX6
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	33019
Uniprot ID	Q7L4S7
NCBI	NP_001171697.1 , NP_001009584.1 , NP_061880.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

