

Product Datasheet

Bovine S100 Calcium Binding Protein A9 (S100A9) ELISA Kit (orb1663697)

Description	Bovine S100 Calcium Binding Protein A9 (S100A9) ELISA Kit
Reactivity	Bovine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.102 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.