

Product Datasheet

OR6K3 antibody (orb166384)

Description	Rabbit polyclonal antibody to OR6K3 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	OR6K3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	37352
Uniprot ID	Q8NGY3
NCBI	NP_001005327.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

