

Product Datasheet

Human RAB1A, Member RAS Oncogene Family (RAB1A) ELISA Kit (orb1663849)

Description Human RAB1A, Member RAS
Oncogene Family (RAB1A) ELISA
Kit

Reactivity Human

**Tested
Applications** ELISA

Note For research use only

Application notes Detect Range: 0.156-10 ng/mL.
Sensitivity: 0.067 ng/mL. Assay
Time: 1-4.5 h. Detection
wavelength: 450 nm. Sample
Volume: 100 µL. Sample Type:
tissue homogenates or other
biological fluids.

Expiration Date Please enquire.