

Product Datasheet

NQ301 (orb16639)

Description	NQ301 is a selective CD45 inhibitor (IC ₅₀ = 200 nM). Exhibits > 200-fold selectivity for CD45 over s...
Note	For research use only
Formula	C ₁₈ H ₁₂ ClNO ₃
MW	325.75
CAS Number	[130089-98-4]
SMILES	CC(=O)C1=CC=C(C=C1)NC2=C(C(=O)C3=CC=CC=C3C2=O)Cl
Expiration Date	12 months from date of receipt.