

Product Datasheet

Dab2 antibody (orb166390)

Description	Rabbit polyclonal antibody to Dab2 (C-term).
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	Dab2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	82312
Uniprot ID	P98078
NCBI	NP_001032994.1 , NP_001008702.1 , NP_075607.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

