

Product Datasheet

IL17RE antibody (orb166393)

Description	Rabbit polyclonal antibody to IL17RE (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	IL17RE
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	74810
Uniprot ID	Q8NFR9
NCBI	NP_001180309.1 , NP_705616.2 , NP_705613.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

