

Product Datasheet

EOMES antibody (orb166395)

Description	Rabbit polyclonal antibody to EOMES (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	EOMES
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	72732
Uniprot ID	O95936
NCBI	NP_001265112.1 , NP_005433.2 , NP_001265111.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

