

Product Datasheet

Bovine Peroxisome Proliferator Activated Receptor Gamma (PPAR γ) ELISA Kit (orb1664000)

Description	Bovine Peroxisome Proliferator Activated Receptor Gamma (PPAR γ) ELISA Kit
Reactivity	Bovine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.115 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.