

Product Datasheet

PTBP1 antibody (orb166401)

Description	Rabbit polyclonal antibody to PTBP1 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PTBP1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	57221
Uniprot ID	P26599
NCBI	NP_002810.1 , NP_114368.1 , NP_114367.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

