

Product Datasheet

Porcine Proinsulin (PI) ELISA Kit (orb1664113)

Description	Porcine Proinsulin (PI) ELISA Kit
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.625-40 pmol/mL. Sensitivity: 0.211 pmol/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.