

## Product Datasheet

### Rat Procollagen I (PCI) ELISA Kit (orb1664249)

|                            |                                                                                                                                                                                                                                                          |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rat Procollagen I (PCI) ELISA Kit                                                                                                                                                                                                                        |
| <b>Reactivity</b>          | Rat                                                                                                                                                                                                                                                      |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                                                                                    |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                    |
| <b>Application notes</b>   | Detect Range: 0.781-50 ng/mL.<br>Sensitivity: 0.36 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm.<br>Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                                                                                          |